**TP BIOLOGIE DE  
SYNTHESE**

3BB 2020

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
200 gel lanes (100 µg) Lot: 1041202 Exp: 2/14  
500 µg/ml Store at -20°C

1.5 ml Gel Loading Dye, Blue (6X) Store at 25°C

**Description:** A number of proprietary plasmids are digested to completion with appropriate restriction enzymes to yield 10 bands suitable for use as molecular weight standards for agarose gel electrophoresis. The digested DNA includes fragments ranging from 0.5–10.0 kilobases (kb). The 3.0 kb fragment has increased intensity to serve as a reference band.

Supplied in: 10 mM Tris-HCl (pH 8.0), 1 mM EDTA.

**Reagents supplied:**  
6X Gel Loading Dye, Blue

**1X Gel Loading Dye, Blue:**  
2.5% Ficoll-400  
11 mM EDTA  
3.3 mM Tris-HCl (pH 8.0@25°C)  
0.017% SDS  
0.015% bromophenol blue

**Preparation:** The double-stranded DNA is digested to completion with appropriate restriction enzymes, phenol extracted and equilibrated to 10 mM Tris-HCl (pH 8.0) and 1 mM EDTA.


**Usage Recommendation:** The 1 kb DNA Ladder was not designed for precise quantification of DNA mass, but can be used for approximating the mass of DNA in comparably intense samples of similar size. The approximate mass of DNA in each of the bands in our 1 kb DNA Ladder is as follows (assuming a 0.5 µg loading):

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**Notes:** All fragments have a 4-base, 5' overhangs that can be end labeled using T4 Polynucleotide Kinase (NEB #M0201) or filled-in using DNA Polymerase I, Klenow Fragment (NEB #M0210) (1). Use α-[<sup>32</sup>P] dATP or α-[<sup>32</sup>P] dTTP for the fill-in reaction.

1 kb DNA Ladder is stable for at least 3 months at 4°C.

For long term storage, store at -20°C. If samples need to be diluted, use TE or other buffer of minimal ionic strength. DNA may denature if diluted in dH<sub>2</sub>O.

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### TSS (Transformation Storage Solution):

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## Q5™ Hot Start High-Fidelity 2X Master Mix

M0494S 002120713071

1-800-432-7799  
info@neb.com  
www.neb.com

# M0494S

100 reactions (50 µl vol) Lot: 0021207

RECOMBINANT Store at -20°C Exp: 7/13

**Description:** The Q5 Hot Start High-Fidelity 2X Master Mix features a high-fidelity, thermostable, hot start DNA polymerase with 3' → 5' exonuclease activity, fused to a processivity-enhancing *Seo7d* domain to support robust DNA amplification. The addition of an aptamer-based inhibitor allows room temperature reaction setup. With an error rate > 50-fold lower than that of *Taq* DNA Polymerase and 6-fold lower than that of *Pyrococcus furiosus* (*Pfu*) DNA Polymerase, Q5 Hot Start High-Fidelity 2X Master Mix is ideal for cloning and can be used for long or difficult amplicons. The convenient master mix formulation is supplied at a 2X concentration. The mix contains dNTPs, Mg<sup>2+</sup> and a proprietary broad-use buffer requiring only the addition of primers and DNA template for robust amplification regardless of GC content. When used at the recommended 1X final concentration, the Q5 Hot Start High-Fidelity Master Mix contains 2 mM MgCl<sub>2</sub>. Q5 Hot Start High-Fidelity 2X Master Mix is unlike typical, lower fidelity PCR enzymes. To determine the optimal annealing temperatures for a given set of primers, use of the **NEB T<sub>m</sub> Calculator** is highly recommended ([www.neb.com/Tmcalculator](http://www.neb.com/Tmcalculator)).

**Please Note:** A precipitate (most noticeable after the first 1–2 freeze/thaw cycles) is not uncommon. To ensure optimal performance, the master mix should be thawed and resuspended prior to use. Stability testing using up to 20 freeze/thaw cycles has shown no negative effect on master mix performance.


**Source:** An *E. coli* strain that carries the Q5 High-Fidelity DNA Polymerase gene.

### Applications:

- High-specificity PCR
- High-fidelity PCR
- Cloning
- Long or difficult amplification
- High-throughput PCR

**Reaction Conditions:** 1X Q5 Hot Start High-Fidelity Master Mix, DNA template and 0.5 µM primers in a total reaction volume of 50 µl.

### Heat Inactivation: No

*Amplification of a variety of human genomic amplicons from *in vivo* high GC content using Q5 Hot Start High-Fidelity 2X Master Mix. All reactions were set up at room temperature, conducted using 30 cycles of amplification and visualized by microfluidic LabChip® analysis.*

### Quality Control Assays

**7 kb Genomic DNA PCR:** 30 cycles of PCR amplification in a 50 µl reaction containing 20 ng genomic DNA with 1X Q5 Hot Start High-Fidelity Master Mix and 0.5 µM of each primer result in the expected 7 kb product.

**20 kb Lambda DNA PCR:** 22 cycles of PCR amplification in a 50 µl reaction containing 10 ng Lambda DNA with 1X Q5 Hot Start High-Fidelity Master Mix and 1.0 µM of each primer result in the expected 20 kb product.

**Hot Start-Specific Genomic DNA PCR:** 25 cycles of PCR amplification in a 25 µl reaction containing 50 ng genomic DNA with 1X Q5 Hot Start High-Fidelity Master Mix in the presence of 0.5 µM of each primer result in the expected 665 bp product, free of non-specific amplification products after pre-incubation at room temperature for 1 hour.

**Note:** Product specifications for individual components in the Q5 Hot Start High-Fidelity 2X Master Mix are available separately.

Oligos : Prom curli sens, Prom curli rev (1µL de primers dilués 1/10<sup>ème</sup>)

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**PCR**

Please note that protocols with Q5 Hot Start High-Fidelity 2X Master Mix may differ from protocols with other polymerases. Conditions recommended below should be used for optimal performance.

**Reaction Setup:**

Q5 Hot Start High-Fidelity 2X Master Mix is inhibited at room temperature, allowing flexible reaction setup (RT or ice).

All components should be mixed prior to use.

COMPONENT	25 µl REACTION	50 µl REACTION	FINAL CONCENTRATION
Q5 Hot Start High-Fidelity 2X Master Mix	1.25 µl	2.5 µl	15 µL 1X
10 µM Forward Primer	1.25 µl	2.5 µl	1,5 µL du stock à 10µM
10 µM Reverse Primer	1.25 µl	2.5 µl	1,5 µL du stock à 10µM
Template DNA	variable	variable	1 µL d'ADNg
Nuclease-Free Water	to 25 µl	to 50 µl	Qsp 30µL

Notes: Gently mix the reaction. Collect all liquid to the bottom of the tube by a quick spin if necessary. Overlay the sample with mineral oil if using a PCR machine without a heated lid.

Transfer PCR tubes to a PCR machine and begin thermocycling.

Q5 Hot Start High-Fidelity 2X Master Mix does not require a separate activation step. Standard Q5 cycling conditions are recommended.

**Thermocycling Conditions for a Routine PCR:**

STEP	TEMP	TIME
Initial Denaturation	98°C	30 seconds
25–35 Cycles	98°C	5–10 seconds
	*50–72°C	10–30 seconds
	72°C	20–30 seconds
Final Extension	72°C	2 minutes
Hold	4–10°C	

\*Use of the NEB T<sub>m</sub> Calculator is highly recommended.

**General Guidelines:**

1. Template:

Use of high quality, purified DNA templates greatly enhances the success of PCR. Recommended amounts of DNA template for a 50 µl reaction are as follows:

DNA	AMOUNT
Genomic	1 ng–1 µg
Plasmid or Viral	1 pg–1 ng

2. Primers:

Oligonucleotide primers are generally 20–40 nucleotides in length and ideally have a GC content of 40–60%. Computer programs such as Primer3 (<http://frodo.wi.mit.edu/primer3>) can be used to design or analyze primers. The best results are typically seen when using each primer at a final concentration of 0.5 µM in the reaction.

3. Mg<sup>2+</sup> and additives:

The Q5 Hot Start High-Fidelity Master Mix contains 2.0 mM Mg<sup>2+</sup> when used at a 1X concentration. This is optimal for most PCR products generated with this master mix.

dNTPs:

The final concentration of dNTPs is 200 µM of each deoxynucleotide in the 1X Q5 Hot Start High-Fidelity Master Mix. Q5 Hot Start High-Fidelity DNA Polymerase cannot incorporate dUTP and is not recommended for use with uracil-containing primers or templates.

5. Q5 Hot Start High-Fidelity DNA Polymerase concentration:

The concentration of Q5 Hot Start High-Fidelity DNA Polymerase in the Q5 Hot Start High-Fidelity 2X Master Mix has been optimized for best results under a wide range of conditions.

6. Denaturation:

Q5 Hot Start High-Fidelity DNA Polymerase does not require a separate activation step.

An initial denaturation of 30 seconds at 98°C is sufficient for most amplicons from pure templates. Longer denaturation times can be used (up to 3 minutes) for templates that require it.

During thermocycling, the denaturation step should be kept to a minimum. Typically, a 5–10 second denaturation at 98°C is recommended for most templates.

7. Annealing:

Optimal annealing temperatures for Q5 Hot Start High-Fidelity 2X Master Mix tend to be higher than for other PCR polymerases. The NEB T<sub>m</sub> Calculator should be used to determine the annealing temperature when using this enzyme. Typically use a 10–30 second annealing step at 3°C above the T<sub>m</sub> of the lower T<sub>m</sub> primer. A temperature gradient can also be used to optimize the annealing temperature for each primer pair.

For high  $T_m$  primer pairs, two-step cycling without a separate annealing step can be used (see note 10).

8. Extension:

The recommended extension temperature is 72°C. Extension times are generally 20–30 seconds per kb for complex, genomic samples, but can be reduced to 10 seconds per kb for simple templates (plasmid, *E. coli*, etc.) or complex templates < 1 kb. Extension time can be increased to 40 seconds per kb for cDNA or long, complex templates, if necessary.

A final extension of 2 minutes at 72°C is recommended.

9. Cycle number:

Generally, 25–35 cycles yield sufficient product.

10. 2-step PCR:

When primers with annealing temperatures  $\geq 72^\circ\text{C}$  are used, a 2-step thermocycling protocol (combining annealing and extension into one step) is possible.

11. Amplification of long products:

When amplifying products > 6 kb, it is often helpful to increase the extension time to 40–50 seconds/kb.

12. PCR product:

The PCR products generated using Q5 Hot Start High-Fidelity 2X Master Mix have blunt ends. If cloning is the next step, then blunt-end cloning is recommended. If T/A-cloning is preferred, the DNA should be purified prior to A-addition, as Q5 Hot Start High-Fidelity DNA Polymerase will degrade any overhangs generated.

Addition of an untemplated -dA can be done with *Taq* DNA Polymerase (NEB #M0267) or Klenow *exo-* (NEB #M0212).

**Companion Products Sold Separately:**

**Q5 Hot Start High-Fidelity DNA Polymerase**

#M0493S	100 units
#M0493L	500 units

**Q5 High-Fidelity DNA Polymerase**

#M0491S	100 units
#M0491L	500 units

**Q5 High-Fidelity 2X Master Mix**

#M0492S	100 reactions
#M0492L	500 reactions

**Q5 Reaction Buffer Pack**

#B9027S	6.0 ml
---------	--------

**Deoxynucleotide Solution Set**

#N0446S	25 $\mu\text{mol}$ of each
---------	----------------------------

**Deoxynucleotide Solution Mix**

#N0447S	8 $\mu\text{mol}$ of each
#N0447L	40 $\mu\text{mol}$ of each

**Magnesium Chloride ( $\text{MgCl}_2$ ) Solution**

#B9021S	6.0 ml
---------	--------

This product is covered by one or more patents.

**Notice to Purchaser:** Nucleic acid-based aptamers for use with thermophilic DNA polymerases are licensed exclusively by New England Biolabs, Inc. from Somalogic, Inc. (See Patent Nos. 5,475,096; 5,670,637; 5,695,249; 5,874,557; and 5,693,502). New England Biolabs, Inc. gives the Buyer/User a non-exclusive license to use the aptamer-based Q5 Hot Start High-Fidelity 2X Master Mix for RESEARCH PURPOSES ONLY. Commercial use of the aptamer-based Q5 Hot Start High-Fidelity 2X Master Mix requires a license from New England Biolabs, Inc. Please contact [busdev@neb.com](mailto:busdev@neb.com) for more information.

This product is licensed from Bio-Rad Laboratories, Inc. under U.S. Pat. Nos. 6,627,424, 7,541,170, 7,760,808, 7,666,645 and corresponding patents in other countries for use only in: (a) standard (non-real time) PCR in the research field only, but not real-time PCR or digital PCR; (b) any *in-vitro* diagnostics application, except for applications using real-time or digital PCR; and (c) any non-PCR applications in DNA sequencing, isothermal amplification and the production of synthetic DNA.

Q5™ is a trademark of New England Biolabs, Inc.

LABCHIP™ is a registered trademark of Galper Life Sciences, Inc., part of PerkinElmer, Inc.


## Clonage SLIC

One-step sequence- and ligation-independent cloning (SLIC):

pIG57 digéré EcoRI (inactiver !)

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### Procedures

1. Digest vector with restriction enzyme(s) overnight, and purify the linearized vector with a commercial PCR purification kit. Elute the DNA with elution buffer or 10 mM TrisCl, pH 8.0-8.5. Do not elute the DNA with water or TE. Measure the concentration of the vector.
2. Amplify your gene of interest by PCR using primers with  $\geq 15$  mer homology extension to the linearized vector end. We usually use 15 bp homology for single fragment cloning, and 20 bp homology for multiple fragment cloning. Purify the linearized vector with a commercial PCR purification kit. Elute the DNA with elution buffer or 10 mM TrisCl, pH 8.0-8.5. Do not elute the DNA with water or TE. Measure the concentration of the insert(s).
3. Mix the linearized vector and insert at a molar ratio of 1:2 in a 1.5 ml tube. An example is shown as follows. (Vector to insert molar ratio of 1:1 to 1:7 works well, but we usually use 1:2 for single fragment cloning, 1:2:2 for multiple fragments cloning. An example of 3 fragments cloning is shown below with vector: insert 1: insert 2 molar ratio is 1:2:2 as shown in Fig. 4A).

	Stock concentration	Volume added	Final concentration
Linearized vector (eg, 5 kb)	100 ng/ $\mu$ l	1 $\mu$ l	10 ng/ $\mu$ l
Insert 1 (PCR product, eg, 1 kb)	40 ng/ $\mu$ l	1 $\mu$ l	4 ng/ $\mu$ l
Insert 2 (PCR product, eg, 1 kb)	40 ng/ $\mu$ l	1 $\mu$ l	4 ng/ $\mu$ l
10X BSA		1 $\mu$ l	1X
10X NEB Buffer 2		1 $\mu$ l	1X
H <sub>2</sub> O		Up to 10 $\mu$ l	

4. Add 0.2  $\mu$ l of T4 DNA polymerase (3 U/ $\mu$ l, NEB) to the mixture and incubate at

room temperature for 2.5 min. We found that 0.5  $\mu\text{l}$  (1.5 U) of T4 DNA polymerase gives the best result, but 0.2  $\mu\text{l}$  gives more than sufficient number of colonies.

5. Put the reaction mixture on ice immediately to stop the reaction and incubate on ice for 10 min.
6. Thaw chemically competent *E. coli* cells on ice for ~10 min. **Utilisez MC4100 *ompR234* « home made »**
7. For single fragment cloning, gently mix the cells with 1-2  $\mu\text{l}$  of the reaction mixture and incubate the cells on ice for 20 min. For multiple fragments cloning, mix the cells with 3-5  $\mu\text{l}$  of the reactant. **Prévoir 3 contrôles : sans ADN, plasmide sans insert, produit de PCR seul**
8. Incubate the cells on ice for 20 min.
9. Heat shock the cells at 42 °C for 45 sec.
10. Incubate the cells on ice for 2 min.
11. Add 900 (or 950)  $\mu\text{l}$  of LB broth to 100 (or 50)  $\mu\text{l}$  of the cells and transfer the cells to 15 ml round-bottom tube.
12. Incubate the cells at 37°C for 1 hr.
13. Plate the cells on agar plates containing suitable antibiotics (eg, 100  $\mu\text{g/ml}$  ampicillin). We usually spread 10 to 20  $\mu\text{l}$  of the cells onto an agar plate to get optimal number of colonies per plate for single fragment cloning, and 100  $\mu\text{l}$  per plate for multiple fragments cloning. In case of 10-20  $\mu\text{l}$  spreading, cells needs to be further diluted with 80-90  $\mu\text{l}$  LB broth for suitable spreading.
14. Incubate the plates at 37°C for 16 hr and analyze the colonies.

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## Miniprep ADN plasmidique

— DNA/RNA  
▲ contaminants

Sample lysis, release of DNA/RNA from cells, tissue, etc.

DNA/RNA is bound to the silica membrane under high-salt conditions  
Interaction between DNA/RNA (hydrate shell is reversibly removed by chaotropic salt) and silica membrane

Contaminants are washed away under high-salt and/or ethanolic conditions to keep the DNA/RNA bound to the membrane

DNA/RNA is eluted in low-salt buffer or water, DNA/RNA is ready to use for downstream applications

## Plasmid DNA purification Protocol-at-a-glance (Rev.08)

1	Cultivate and harvest bacterial cells	 	11,000 × g, 30 s
2	Cell lysis		250 µL Buffer A1 250 µL Buffer A2 RT, up to 5 min 300 µL Buffer A3
3	Clarification of the lysate	 	11,000 × g, 5–10 min
4	Bind DNA	 	Load supernatant 11,000 × g, 1 min
5	Wash silica membrane	 	(Optional: 500 µL Buffer AW: RT or 50 °C) 600 µL Buffer A4 11,000 × g, 1 min
6	Dry silica membrane	 	11,000 × g, 2 min
7	Elute DNA	 	50 µL Buffer AE RT, 1 min 11,000 × g, 1 min

## 5 NucleoSpin® Plasmid/Plasmid (NoLid) protocols

### 5.1 Isolation of high-copy plasmid DNA from *E. coli*

#### Before starting the preparation:

- Check if Wash Buffer A4 was prepared according to section 3.

#### 1 Cultivate and harvest bacterial cells

Use 1–5 mL of a saturated *E. coli* LB culture, pellet cells in a standard benchtop microcentrifuge for **30 s** at **11,000 x g**. Discard the supernatant and remove as much of the liquid as possible.

**11,000 x g,  
30 s**

*Note: For isolation of low-copy plasmids refer to section 5.2.*

#### 2 Cell lysis

Add **250 µL Buffer A1**. Resuspend the cell pellet completely by vortexing or pipetting up and down. Make sure no cell clumps remain before addition of Buffer A2!

**+ 250 µL A1**

**Resuspend**

*Attention: Check Buffer A2 for precipitated SDS prior to use. If a white precipitate is visible, warm the buffer for several minutes at 30–40 °C until precipitate is dissolved completely. Cool buffer down to room temperature (18–25 °C).*

**+ 250 µL A2**

**Mix**

**RT, 5 min**

Add **250 µL Buffer A2**. Mix gently by inverting the tube **6–8 times**. Do not vortex to avoid shearing of genomic DNA. Incubate at **room temperature** for up to 5 min or until lysate appears clear.

Add **300 µL Buffer A3**. Mix thoroughly by inverting the tube **6–8 times**. Do not vortex to avoid shearing of genomic DNA!

**+ 300 µL A3**

**Mix**

#### 3 Clarification of lysate

Centrifuge for **5 min** at **11,000 x g** at room temperature.

Repeat this step in case the supernatant is not clear!

**11,000 x g,  
5–10 min**

#### 4 Bind DNA

Place a NucleoSpin® Plasmid/Plasmid (NoLid) Column in a Collection Tube (2 mL) and decant the supernatant from step 3 or pipette a maximum of 750 µL of the supernatant onto the column. Centrifuge for **1 min** at **11,000 x g**. Discard flow-through and place the NucleoSpin® Plasmid/Plasmid (NoLid) Column back into the collection tube.

**Load  
supernatant**

**11,000 x g,  
1 min**

Repeat this step to load the remaining lysate.

#### 5 Wash silica membrane

*Recommended: If plasmid DNA is prepared from host strains containing high levels of nucleases (e.g., HB101 or strains of the JM series), it is strongly recommended performing an additional washing step with 500 µL Buffer AW, optionally preheated to 50 °C, and centrifuge for 1 min at 11,000 x g before proceeding with Buffer A4. Additional washing with Buffer AW will also increase the reading length of DNA sequencing reactions and improve the performance of critical enzymatic reactions.*

**Optional:  
+ 500 µL AW**

**11,000 x g,  
1 min**

Add **600 µL Buffer A4** (supplemented with ethanol, see section 3). Centrifuge for **1 min** at **11,000 x g**. Discard flow-through and place the NucleoSpin® Plasmid/Plasmid (NoLid) Column back into the empty collection tube.

**+ 600 µL A4**

**11,000 x g,  
1 min**

#### 6 Dry silica membrane

Centrifuge for **2 min** at **11,000 x g** and discard the collection tube.

*Note: Residual ethanolic wash buffer might inhibit enzymatic reactions.*

**11,000 x g,  
2 min**

#### 7 Elute DNA

Place the NucleoSpin® Plasmid/Plasmid (NoLid) Column in a 1.5 mL microcentrifuge tube (not provided) and add **50 µL Buffer AE**. Incubate for **1 min** at **room temperature**. Centrifuge for **1 min** at **11,000 x g**.

**+ 50 µL AE**

**RT, 1 min**

*Note: For more efficient elution procedures and alternative elution buffer (e.g., TE buffer or water) see section 2.4.*

**11,000 x g,  
1 min**

# Annexes

## Sequence Promoteur curli P<sub>curli</sub> csgBA

<http://www.ecogene.org/?q=ecodownload/sequence>

<http://www.ecogene.org/?q=gene/EG12621>

>E. coli K-12 MG1655 U00096.3 (1103197 to 1103950 = 754 bp)

```
GATGAAACCCCGCTTTTTTATTGATCGCACACCTGACAGCTGCCTCTAAAATAGAAGCACCAGAAGTACTGACA
GATGTTGCACTGCTGTGTGTAGTAATAAATCAGCCCTAAATGGGTAAAATATAAACTAATGGATTACATCTGAT
TTCAATCTAGCCATTACAAATCTTAAATCAAGTGTTAAACATGTAACCTAAATGTAACCTCGTTATATTTAAATGTT
AACCTTAAGGTTTTATTAAGTTTAGAAATGATAGAAAAGTTGTACATTTGGTTTTTATTGCACAATTTAAAAAA
TCATACAAATGGTGATAACTACTAATAATGCATATAAAAAATATTTGGGTGTAGTCCTTTGTGCATGTAAAACG
TTCTTGTTTTTTCTCCACACCTCCGTGGACAATTTTTTACTGCAAAAAGACGAGGTTTGTACGGCTTGTGCGCA
AGACATATCGCAGCAATCAGCGACGGGCAAGAAGAATGACTGTCTGGTGTCTTTTGATAGCGGAAAACGGAGATT
TAAAAGAAAACAAAATATTTTTTTCGTTAGATAACAGCGTATTTACGTGGGTTTTAATACTTTGGTATGAACTAA
AAAAGAAAATACAACGCGCGGGTGTAGTTATTTAAAAATATTTCCGCAGACATACTTTCCATCGTAACGCAGCGTT
AACAAAATACAGTTTCGTTAACACCAAGTTGAAATGATTTAATTTCTTAAATGTACGACCAGGTCCAGGGTGA
CAAC
```

## Carte pSB1K3

<http://parts.igem.org/Part:pSB1K3>

<http://parts.igem.org/cgi/partsdb/puttext.cgi>

## Registry of Standard Biological Parts

main page
design
experience
information
part tools
edit

### Part:pSB1K3

Designed by: Austin Che Group: iGEM07\_Example (2008-08-08)


Not Released
Sample Not in stock
★ 1 Registry Star
Not Used
Get This Part

#### High copy BioBrick assembly plasmid

This is a high copy plasmid for use in assembly of BioBrick® standard biological parts.

Sequence and Features

Subparts | [Ruler](#) | [SS](#) | [DS](#)      Length: 2204 bp      [View plasmid](#) | [Get part sequence](#)

Assembly Compatibility: 10 12 21 23 25 1000

Parameters


chassis	
copies	100-300
insert	None
mc	BioBrick
origin	pMB1
resistance	K

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